

A New Algorithm for FAIMS Data Analysis with Accurate In-depth Quantitative Profiling

Zia Rahman; Yandong Zhu; Zheng Chen; Jonathan R. Krieger. Baozhen Shan
Bioinformatics Solutions Inc., Waterloo, ON

EXPANDED ABSTRACT


Peptide identification by data-dependant mass spectrometry relies on the sequential selection of precursor ions as they elute from the chromatographic column. Complex samples represent a challenge due to the presence and occurrence of chimeric tandem mass spectra from co-eluting ions. Gas-phase ion fractionation opens new perspectives to overcome sample complexity. High-field asymmetric waveform ion mobility spectrometry (FAIMS) provides an important resource to combat this issue, however, quantitative analysis of FAIMS data remains a challenge. Here we present a new algorithm, including 4-dimensional feature detection and alignment for FAIMS data analysis. Preliminary testing indicates that algorithm enabled deep proteome coverage with high accuracy.

The algorithm contains four components:

- FAIMS holds the promise of delivering gas-phase fractionation by stepping through multiple compensation voltages (CVs), deconvoluting overlapping peptide signals, and removing contaminants. Thus, peptide precursor features are detected by grouping each CV from the run together allowing for higher precursor resolution
- Based on precursor m/z , retention time, and CV, multiple features are associated with each MS/MS spectrum. This association further deconvolutes overlapping precursors for in-depth identification
- Label free quantification is achieved by performing 4D-feature alignment across different LC-MS/MS runs. Peptide intensity is calculated as summation of values in all CVs.
- TMT quantification is achieved by excluding chimeric spectra, and by also separating MS/MS spectra from distinct CV values resulting in higher accuracy quantification.

The algorithm was implemented in PEAKS XPro software platform and tested with multiple published available data sets. With one dataset (PXD009547), the new algorithm suggested that FAIMS aids in reducing the percentage of chimeric spectra from 39% to 25%. This reduction in chimerism leads to significant benefits in improving proteome coverage. A 36% gain in the number of protein identifications (2 unique peptides) compared with traditional LC-MS/MS analyses for the same number of injections. The reduction in complexity also has improving effects on quantification. With a second dataset, an increase of 1100 protein groups were quantified using FAIMS, with noticeable increases in feature detection and alignment correlations, and decreased in standard deviation across injection replicates. Together, this suggests that the new PEAKS algorithm provides accurate and sensitive quantitative analysis with FAIMS data.


INTRODUCTION



Hebert, A. S., Prasad, S., Belford, M. W., Bailey, D. J., McAlister, G. C., Abbatiello, S. E., et al. (2018). Comprehensive Single-Shot Proteomics with FAIMS on a Hybrid Orbitrap Mass Spectrometer. *Analytical Chemistry*.

Principles of FAIMS

- FAIMS (high-Field Asymmetric waveform Ion Mobility Spectrometry) allows for the separation of ions based on their mobility in gas phase
- Ions get trapped by a carrier gas between two parallel electrodes. Adding a high frequency waveform (DV) to the inner electrode while the outer electrode remains grounded causes ions to be trapped




- A compensation voltage (CV) is added to the electrode to counteract the DV allowing ions to be transmitted through the FAIMS device into the mass spectrometer (MS)
- Different CVs will cause different ions to move through the FAIMS device into the MS.
- Ions of different charge states will distribute different as a function of CV
- This represents a second separation step to simplify the peptide mixture entering the mass spectrometer to improve detection

Cooper, H. J. (2016). To What Extent is FAIMS Beneficial in the Analysis of Proteins? *Journal of the American Society for Mass Spectrometry*.

Purves, R. W., & Guevremont, R. (1999). Electrospray Ionization High-Field Asymmetric Waveform Ion Mobility Spectrometry–Mass Spectrometry. *ACS Publications*

Pfammatter, S., Bonneil, E., McManus, F. P., Prasad, S., Bailey, D. J., Belford, M., et al. (2018). A Novel Differential Ion Mobility Device Expands the Depth of Proteome Coverage and the Sensitivity of Multiplex Proteomic Measurements. *Molecular & Cellular Proteomics : MCP*, 17(10), 2051–2067.

PEAKS SOLUTIONS FOR 4D FAIMS-LC-MS/MS DATA



- PEAKS X PRO Studio and PEAKS Online are complete proteomics software solutions that now includes full support for FAIMS 4D Proteomics Projects, including:
 - FAIMS DDA
 - FAIMS DIA*
 - FAIMS SILAC
 - FAIMS TMT
 - Multiple CV values in a single run



*direct DIA in Studio, Library search in both Studio and Online/

Both versions of PEAKS employ a new algorithm specifically geared towards FAIMS data.


FAIMS Algorithm Feature	Application/Function
Individual CVs are grouped together	higher precursor resolution and better feature detection
Multiple overlapping feature associations per MS/MS spectrum	increased in-depth protein identification
Performs 4D feature alignment across different LC-MS/MS Runs	peptide intensity is calculated as a summation of values in all CVs for better and more accurate quantification
Specific to TMT, chimeric spectra are excluded and MS/MS spectra from distinct CV values are separated	higher quantification accuracy

IMPROVEMENTS IN ANALYSIS OF FAIMS vs. NON-FAIMS DATA


Using data from Pfammatter S. *et al*¹, which specifically compared non-FAIMS to FAIMS data and used PEAKS 8.5 for analysis, we sought to confirm and show how improvements in the algorithm have led to a marked increase in protein identifications



Protein Identifications (>= 2 Unique Peptides)



Unique Peptide Sequences




Re-analyzing this data set with the new PEAKS X PRO shows substantial gains in protein identifications for data generated both with and without FAIMS, but specifically the largest gains in sensitivity and identifications of both peptides and proteins comes with data generated using FAIMS.

¹Pfammatter, S., Bonneil, E., McManus, F. P., Prasad, S., Bailey, D. J., Belford, M., et al. (2018). A Novel Differential Ion Mobility Device Expands the Depth of Proteome Coverage and the Sensitivity of Multiplex Proteomic Measurements. *Molecular & Cellular Proteomics : MCP*, 17(10), 2051–2067.

PEAKS STUDIO X PRO and PEAKS ONLINE offer Consistent Results


In order to assess whether there are differences in the identification of peptides and proteins between PEAKS Studio X Pro and PEAKS Online (v.1.4), we carried out a direct comparison using the FAIMS triplicate injections of HEK293T cells with multiple CV values run in DDA mode from the same dataset.





As illustrated in the Venn Diagrams, we observe >95% overlap between PEAKS Studio X Pro and PEAKS Online with DDA data.

CV and LFQ ANALYSIS

Further examining the same dataset, proteins identified in all replicate injections were quantified. CV values remain low both with and without FAIMS, across 3 replicates in each condition. However, quantification across the FAIMS data shows a broader distribution, likely due to several factors including more proteins quantified; more accurate quantification; the reduction in chimeric spectra while using FAIMS, and improvement of sensitivity including minimizing ion suppression. With the same dataset between PEAKS Studio X Pro and PEAKS Online, we obtain very similar results of label free quantification



TMT, Phosphorylation, and FAIMS with PEAKS



- With the clear beneficial advantages of FAIMS in reducing sample complexity and increasing identifications, its use and presence in the field will become more pronounced. Furthermore, due the increased throughput of isobaric labelling by Tandem Mass Tags, we have introduced the ability to do TMT quantification with the PEAKS line of products.
- In order to assess the ability of the PEAKS tools to identify and quantify phosphopeptides with isobaric tags, we took a published dataset from Schweppe *et al*, and compared the numbers of unique peptide identifications across a number of conditions (Synchronous Precursor Selection MS3 (SPS-MS3) without FAIMS, High Resolution MS2 (HR-MS2) without FAIMS, or HR-MS2 with various FAIMS CVs as indicated.
- We observe results that benchmark extremely well against this published dataset.
- Isobaric labels can then be used for accurate quantification

Schweppe, D. K., Rusin, S. F., Gygi, S. P., & Paulo, J. A. (2019). Optimized Workflow for Multiplexed Phosphorylation Analysis of TMT-Labeled Peptides Using High-Field Asymmetric Waveform Ion Mobility Spectrometry. *Journal of Proteome Research*, 19(1), 554–560.

Summary

- PEAKS Studio X Pro and PEAKS Online v1.4 can represent a new software solution for analyzing FAIMS Pro Data
- All types of FAIMS data can be analyzed, including:
 - DDA, DIA, TMT, SILAC
- Multiple CVs in same files can be used directly with PEAKS without the need for conversion or splitting data files
 - This also means that individual CV values do not need to be run independently of each other. CV switching on the FAIMS Pro device is relatively quick, and thus a combination approach has been shown to be the best
- PEAKS Studio X Pro can directly identify peptides and proteins from DIA FAIMS experiments if a library is not available. This is done accurately using our well known *de novo* sequencing algorithm
- For a demonstration or a free trial, feel free to visit us or contact us (peaks@bioinfor.com)